

Product Comparison



	Model M	Model A
FDA-cleared	✓	No
Partially reimbursable by Medicare and private insurance	Yes, for patients who qualify	No
Armrests	✓	No
iPhone app	No	✓
Distance with batteries fully charged	15 miles	12 miles
Maximum speed	5.5mph	5.5mph
Back support options	Rehab low, med, tall, deep	low, mid foam filled
Seat cushion	Rehab (skin protecting)	General use (foam)
Standard back post expansions	✓	No
For transporting in a vehicle	4-point tie-down	No
Pressure relief	✓	No

Specifications for both models	Model M	Model A
5-year product warranty	✓	✓
4WD	✓	✓
3" Obstacle clearance Over curbs, rocks; through snow, rain	✓	✓
10 degrees of incline/decline	✓	✓
Power seat slide - forward 5.9"	✓	✓
Both arms rotate up to get in and out	✓	✓
Mouse controller - multiple types	✓	✓
Tight turning radius of 28"	✓	✓
Electro-magnetic braking system	✓	✓
220lb max weight limit	✓	✓
Attach accessories such as a cupholder, iPhone, iPad, flashlight	✓	✓
List price	\$12,499	\$12,999
Approximate cost with insurance	\$7,000	\$12,999
Sold through various health providers around the country	Yes	Yes
Sold online	No	Yes (spinlife.com) Check website for pricing

Note: Caution Model M is FDA-cleared and Federal law restricts this device for sale by or on the order of a physician. The Model A is not considered a medical device and has not been submitted to the FDA for review or clearance.

740-01951-MKT-02