


WHILL Model M Fitting, Customization, and Specifications

With a wide degree of adjustability and ergonomically designed seating, as well as broad compatibility with accessories, the Model M can be customized and accessorized to meet the needs of a wide range of users.



Controller Options


The ergonomic controller comes in multiple shapes, making it appropriate for a variety of dexterity levels.



MOUSE



EASY GRIP



Right or left hand side

The directional controller can be configured on the right or left hand side.

Back Support Options



LOW



MID



TALL



DEEP

Pelvic Stabilizing Belt

Available in medium and large



Performance Specifications

Driving Range (full battery charge)	15 miles
Maximum Speed	5.5 MPH
Speed Settings	High, Mid, Low
Drive System	Two Motor Drive System
Incline Capability	10 Degrees
Seat Slide Range (electric)	5.9"

Controller	WHILL Control System
Turning Radius	28"
Max Obstacle Height	3"
Chair Ground Clearance	3.5"
Braking System	Electromagnetic
Max Weight Capacity	220 lb

Dimensional Specifications

Chair Weight	250 lb
Chair Width	23.6"
Chair Length	35" - 39"

Front Wheel Diameter	9.8"
Rear Wheel Diameter	12.5"
Rear Tires	Solid

Other Specifications

Charging Time	8 - 9 Hours
Battery	12V 50Ah (per battery; 2 total batteries)
Charger	6A Charger

Operating Conditions	-13°F to 122°F
Storage Conditions	-40°F to 149°F
Tail Lamps (2)	Red Lights

WHILL, Inc. 285 Old County Road Suite 6 San Carlos, CA 94070 (844) 699 - 4455

CAUTION: Federal law restricts this device to sale by or on the order of a physician
740-01529-MKT-02 3/2017

www.whill.us